



MOUSE ANTI ADAMTS13 MONOCLONAL ANTIBODY

Catalog No.	LNM-KR-051
Clone	A10
Concentration	1 mg/ml
Volume	100 UL
Subclass	IgG2bk
Immunogen	recombinant ADAMTS13
Specificity	ADAMTS13 [disintegrin-like domain]
Applications	Suitable for immunoblotting (non-reducing condition), Immunohistochemistry, ELISA (capture antibody)
Format	Purified immunoglobulin from BALB/C Mouse Ascites Ammonium sulfate fractionation Affinity chromatography
Formulation	Liquid in 15mM Phosphate buffer containing 0.15M NaCl, 0.09% sodium azide, pH 7.2
Storage/Handling	2 to 8°C, Maintain at -70°C for up to 12 months. Avoid repeated freeze/thaw cycle.
Reference	Uemura M et al., Blood., 106, 922-924 (2005)

For research use only; not for use as a diagnostic



抗 ADAMTS13 モノクローナル抗体

品 番	KR-051
クローン	A10
型	IgG2b, κ
免疫抗原	recombinant ADAMTSA13
特異性	ADAMTS13 (Disintegrin-like domain)
由 来	BALB/Cマウス腹水
精製法	硫酸塩析, r-プロテインA
溶 液	PBS (0.09%アジ化ナトリウム含む)
数 量	100 μ g
保 存 法	4~8℃ 長期保存の場合は、-70℃以下にて凍結保存して下さい。 凍結融解は出来るだけ避けて下さい。
備 考	WB (非還元条件)、IH、ELISA (固相抗体)
Reference	Uemura M et al., Blood., 106., 922-924.(2005)



コスモ・バイオ株式会社