




ALUM

Form	Aluminium hydroxide hydrate gel suspension
Volume	100mg : 0.9% NaCl solution (20mg/ml)
Molecular Formula	$\text{Al}(\text{OH})_3 \cdot x\text{H}_2\text{O}$ gel suspension
Molecular Weight	78.00 (anhydrous basis)
Immunity grade	adjuvant use
Storage	Store below 4°C (below -20°C for prolonged storage)



Allergic antigen

LG-0009	DNP-Ascaris
LG-0533	DNP-Mite-Df
LG-0544	DNP-Mite-Dp
LG-0528	DNP-Cedar Pollen-Cj
LG-0017	DNP-BSA
LG-3017	DNP-BSA (30)
LG-0024	DNP-OA (Ovalbumin)
LG-0047	DNP-Casein
LG-0067	DNP-LG (BETA-Lactoglobulin)
LG-0089	DNP-KLH
LG-0317	BPO-BSA
LG-0577	BPO-BGG
LG-1109	TNP-Ascaris
LG-1117	TNP-BSA
LG-5009	Ascaris Extract
LG-5339	Mite Extract-Df
LG-5449	Mite Extract-Dp
LG-5559	Mite Extract-Tp
LG-5229	Cedar Pollen Extract
LG-5280	Cedar Pollen Extract-Cj
LG-5779	Hinoki Cypress, Pollen Crude Extract
LG-5889	Ragweed Pollen Extract

Adjuvant

LG-6000	ALUM
---------	------

For research use only. Not for clinical diagnosis.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619

ALUM (水酸化アルミニウムゲル)

I. 内容

lot No. 745022

性 状	Aluminum hydroxide hydrate $[Al(OH)_3 \cdot xH_2O]$ gel suspension
規 格	100mg : 0.9% NaCl 溶液に懸濁 (20mg/ml) 免疫グレード : アジュバンド用

II. 保存上の注意

4℃以下 (長期間の場合は-20℃以下) で保存して下さい。

製造元

総発売元

 株式会社 エル・エス・エル

コスモ・バイオ株式会社

新製品案内

- ☆Anti (MAP) Cry j-1
カタログNo. LB-5201 (LSL)
- ☆Anti (MAP) Cry j-2
カタログNo. LB-5202 (LSL)
- ☆Ragweed Pollen Extract
カタログNo. LG-5889 (LSL)
- ★Cedar Pollen Extract-Cj
カタログNo. LG-5280 (LSL)
- ★FTPP (フォルマリン処理百日咳死菌) 懸濁液
カタログNo. LG-7000 (LSL)