



COSMO BIO CO., LTD.
Inspiration for Life Science

Monoclonal Antibody

to Human Chorionic Gonadotropin(MCG-033)

CATALOG NUMBER:	NBT-MCG-033
QUANTITY:	1mg
CLONE:	CG-033
SOURCE OF ANTIGEN:	Human Chorionic Gonadotropin
SOURCE:	BALB/C Mouse implantation ascites
FORMAT:	Purified Protein
CONCENTRATION:	10mg/ml
ISOTYPE:	IgG1k
PURIFY:	Protein A
PURITY:	1 band (IgG1k) Cellulose acetate Electrophoresis
PRESENTATION:	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles

For research use only; not for use as a diagnostic.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp