



COSMO BIO CO., LTD.
Inspiration for Life Science

Monoclonal Antibody to D Dimer (DF-041)

CATALOG NUMBER:	NBT-MDF-041
LOT NUMBER:	010411-1
QUANTITY:	1mg
CLONE:	DF-041
SUBCLASS:	IgG1k
SOURCE:	BALB/C Mouse implantation ascites
FORMAT:	Purified Protein
CONCENTRATION:	10mg/ml
PURIFY:	Protein A
PURITY:	1 band (IgG1k) Cellulose acetate Electrophoresis
PRESENTATION:	Liquid in 0.15M NaCl 10mM PB, pH7.2, 0.1% Sodium Azide
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.


For research use only; not for use as a diagnostic.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp



Nippon Bio-Test Laboratories Inc.

株式会社日本バイオテスト研究所
〒115-0022 東京都国分寺市東戸倉1-15-3 直(042)321-0921

抗ヒトDDダイマーモノクローナル抗体

カタログ No.	#MDF-041	#MDF-042	#MDF-043	#MDF-044
クローン	DF-041	DF-042	DF-043	DF-044
サブクラス	IgG _{1κ}	IgG _{1κ}	IgG _{2ακ}	IgG _{1κ}
Lot No.	010411-1	010418-1	010412-1	010413-1

- 抗体量 各1 mg
- 蛋白濃度 10mg/ml (OD280 1%=14.0)
- 由来 BALB/Cマウス移植腹水
- 精製法 プロテインA, 硫酸塩析
- 溶液 0.15M NaClを含む10mM PB, pH7.2
- 精製純度 セルロースアセテート膜電気泳動にて単一な抗体バンドを示す。
- 防腐剤 0.1%アジ化ナトリウム
- 保存法 4~8°C
長期保存の場合は-20°C以下にて凍結保存が有効です。
凍結融解はできるだけ避けて下さい。
- 備考

株式会社 日本バイオテスト研究所
東京都国分寺市東戸倉1丁目15-3
TEL 042 321-0921/FAX 042 323-7025

コスモ・バイオ株式会社