



Monoclonal Antibody to Human α Fetoprotein (PG-012)

CATALOG NUMBER:	NBT-PG-012
QUANTITY:	5ml (0.5mg/ml)
CLONE:	PG-012
SOURCE OF ANTIGEN:	Human α Fetoprotein
SOURCE:	Rabbit Searum
FORMAT:	Purified Protein
CONCENTRATION:	6.3mg/ml
ISOTYPE:	IgG
PURIFY:	Protein A
PURITY:	1 band Cellulose acetate Electrophoresis
PRESENTATION:	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles

For research use only; not for use as a diagnostic.



NBT#PA012

品番. PG-012

試験成績書

製品名: 抗ヒトAFPウサギIgG 分画 Lot No. _____

試験項目	試験結果
①特異抗体量	0.5 mg/ml (Becker 法)
②蛋白濃度	mg/ml (吸光度法)
③容 量	5.0 ml
④免疫動物	ウサギ・ヤギ・ヒツジ・ニワトリ
⑤純度検定 a: 免疫電気泳動法 b: オクタロニー法 (多検体及び希釈試験)	AFP以外の沈降線は観察されない 不純抗体は認められない
⑥貯 法	2~8℃
⑦性 状	0.15M NaCl, 0.09%NaN ₃ を含む 10mMリン酸緩衝液 (pH7.2)に溶解
⑧備 考	

東京都国分寺市東戸倉1-15-3
株式会社日本バイオテスト研究所
T e l (042) 321-0921
F a x (042) 323-7025

コスモ・バイオ株式会社