



Monoclonal Antibody to C Reactive Protein (PG-022)

CATALOG NUMBER:	NBT-PG-022
QUANTITY:	5ml (0.5mg/ml)
CLONE:	PG-022
SOURCE OF ANTIGEN:	C Reactive Protein
SOURCE:	Rabbit Searum
FORMAT:	Purified Protein
CONCENTRATION:	9.9mg/ml
ISOTYPE:	IgG
PURIFY:	Protein A
PURITY:	1 band Cellulose acetate Electrophoresis
PRESENTATION:	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles

For research use only; not for use as a diagnostic.



品番. PG-022

試験成績書

製品名：抗CRPウサギIgG分画 Lot.No. _____

試験項目	試験結果
①特異抗体量	2.0 mg/ml (Becker 法)
②蛋白濃度	mg/ml (吸光度法)
③容 量	5.0 ml
④免疫動物	ウサギ・ヤギ・ヒツジ・ニワトリ
⑤純度検定 a : 免疫電気泳動法	CRP以外の沈降線は観察されない
⑥貯 法	2 ~ 8 °C
⑦性 状	0.15M NaCl, 0.09%NaN ₃ を含む 10mMリン酸緩衝液(pH 7.2)に溶解

東京都国分寺市東戸倉1-15-3
株式会社日本バイオテスト研究所
T e l (0 4 2) 3 2 1 - 0 9 2 1
F a x (0 4 2) 3 2 3 - 7 0 2 5

コスモ・バイオ株式会社