



COSMO BIO CO., LTD.
Inspiration for Life Science

Monoclonal Antibody to Human Hemoglobin (PS-041)

CATALOG NUMBER:	NBT-PS-041
QUANTITY:	10ml
CLONE:	PS-041
SOURCE OF ANTIGEN:	Human Hemoglobin
SOURCE:	Goat
FORMAT:	Searum
CONCENTRATION:	98.0 mg/ml
PRESENTATION:	Liquid in 10mM phosphate buffer with 0.15M NaCl , pH7.4, 0.1% Sodium Azide
STORAGE/HANDLING:	Maintain at 4 °C-8 °C (-20 °C for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles

For research use only; not for use as a diagnostic.



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO EKIMAE BLDG, 2-20, TOYO 2CHOME, KOTO-KU, TOKYO 135-0016, JAPAN
TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 e-mail: export@cosmobio.co.jp

www.cosmobio.co.jp

