



Anti-Mouse KIAA0202 Polyclonal Antibody, Rabbit

Mouse KIAA0202 (mKIAA0202) protein is a homologue of human KIAA0202 (ref 1, 2). KIAA0202 is identical to Septin 8 (ref. 3). Rabbit anti mKIAA0202 antibody is raised against GST-fused recombinant protein (GX0010) containing following sequence.

(192 amino acids)

AVVGSTEEVKVGNKLVRRARQYPWGVVQVENENHCDFVKLREMLIRVNMEDLREQTHSRHYELY
RRCKLEEMGFQDSDGDSQPFSLQETYEAKRKEFLSELQRKEEEMRQMFVNKVKETELELKEKE
RELHEKFEHLKRIHQEEKRKVEEKRRLEEEETNAFNCRKAAMEALQSQUALHATSQQPLRKDKDK
KN

Use in any such application is responsibility of the user.

Preparation: Purified by the affinity chromatography from rabbit antisera and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0010. This antibody detects endogenous levels of mKIAA0202 protein in several cell types. It also recognizes human KIAA0202 protein. Details are shown in the InGaP database. (<http://webcreate.kazusa.or.jp/create/index.jsp>)

Application: Western blotting (1 : 1,000), Immunohistochemistry (frozen sections, paraffin-embedded sections), Immunoprecipitation.

Storage: Keep at -20°C

Package: 50 μg (250 μL)

Product Code: MK02020505

References:

- 1) Hara, Y. *et al.*: *DNA Res.*, 10, 129 (2003).
- 2) Koga, H. *et al.*: *DNA Res.*, 11, 293 (2004).
- 3) Nagata, K. *et al.*: *J. Biol. Chem.*, 279, 55895 (2004)

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619