



Anti-Mouse KIAA0696 Polyclonal Antibody, Rabbit

Mouse KIAA0696 (mKIAA0696) protein is a homologue of human KIAA0696 (ref 1, 2). KIAA0696 is identical to F-box and WD-40 domain protein 11 (ref. 3). Rabbit anti mKIAA0696 antibody is raised against GST-fused recombinant protein (GX0195) containing following sequence.

(108 amino acids)

CLRVLEGHEELVRCIRFDNKRIVSGAYDGKIKVWDLQAALDPRAPASTLCLRTLVEHSGRVFRLQ
FDEFQIISSSHDDTILIWDFLNVPPSAQNETRSPSRITYTYISR

Use in any such application is responsibility of the user.

Preparation: Purified by Protein A affinity chromatography from rabbit antisera and solubilized in PBS containing with 50% glycerol and 0.02% of NaN₃.

Specificity: Specific to recombinant protein GX0195. This antibody detects endogenous levels of mKIAA0696 protein in cerebellar purkinje and molecular layer cells.

Application: Western blotting (1 : 1,000), Immunohistochemistry

Storage: Keep at -20°C

Package: 100 µg (200 µL)

Product Code: MK06960310

References:

- 1) Hara, Y. *et al.*: *DNA Res.*, 10, 129 (2003).
- 2) Koga, H. *et al.*: *DNA Res.*, 11, 293 (2004).
- 3) Cenciarelli, C. *et al.*: *Curr Biol.*, 9, 1177 (1999).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619