



Anti-Mouse KIAA1201 Polyclonal Antibody, Rabbit

Mouse KIAA1201 (mKIAA1201) protein is a homologue of human KIAA1201 (ref 1, 2). KIAA1201 is identical to GRAM domain containing 1B. Rabbit anti mKIAA1201 antibody is raised against GST-fused recombinant protein (GX0439) containing following sequence.

(370 amino acids)

```
CEEIPIEENEVNDSSSKSSIETKPDASPQLPKKSITNSTLTSTGSSEAPVSFDGLPLEEEVMEGDG  
SLEKELAIDNIGEKIEIMAPVTSPSLDFNDNEDIPTELSDSSDTHDEGEVQAFYEDLSGRQYVNE  
VFNFSVDKLYDLLFTNSPFLRDFMEQRRFSDIIFHPWKKEENGNQSRVILYTITLTNPLAPKTATVR  
ETQTMKASQESECYVIDAEVLTHDVPYHDYFYTINRYTLTRVARNKSRLRVSTELRYRKQPWG  
FVKTFIEKNFWSGLEDYFRHLETETKTESTYLAEIHRQSPKEKASKSSAVRRRKRPHAHLRVPH  
LEEVMSPVTTPTDEDVGHRIKHVAGSTQTRHIPEDTPDGFHL
```

Use in any such application is responsibility of the user.

Preparation: Purified by Protein A affinity chromatography from rabbit antisera and solubilized in PBS containing with 50% glycerol and 0.02% of NaN₃.

Specificity: Specific to recombinant protein GX0439. This antibody detects endogenous levels of mKIAA1201 protein in cerebellar interneuron cells.

Application: Western blotting (1 : 1,000), Immunohistochemistry

Storage: Keep at -20°C

Package: 100 µg (200 µL)

Product Code: MK12010310

References:

- 1) Hara, Y. *et al.*: *DNA Res.*, 10, 129 (2003).
- 2) Koga, H. *et al.*: *DNA Res.*, 11, 293 (2004).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO Co., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619