



Anti-Mouse KIAA1465 Polyclonal Antibody, Rabbit

Mouse KIAA1465 (m KIAA1465) protein is a homologue of human KIAA1465 (ref 1, 2). KIAA1465 is identical to immunoglobulin superfamily containing leucine-rich repeat 2. Rabbit anti mKIAA1465 antibody is raised against GST-fused recombinant protein (GX0968) containing following sequence.

(145 amino acids)

LVLATVPLLGAACCHLLAKHPGKPYRLILRPQAPDPMEKRIAADFDPRASYLESEKSYPARGEAG
GEEPEEVPEEGLDEDVEQGDPSGDLQREESLAGCSLVESQSKANQEEFEAGSEYSDDLPLGAE
AVNIAQEINGNYRQTAG

Use in any such application is responsibility of the user.

Preparation: Purified by Protein A affinity chromatography from rabbit antisera and solubilized in PBS containing with 50% glycerol and 0.02% of NaN_3 .

Specificity: Specific to recombinant protein GX0968. This antibody detects endogenous levels of mKIAA1465 protein in several cell types. Details are shown in the InGaP database.

(<http://webcreate.kazusa.or.jp/create/index.jsp>)

Application: Western blotting (1 : 1,000), Immunoprecipitation, Immunohistochemistry.

Storage: Keep at -20°C

Package: 100 μg (200 μL)

Product Code: MK14650310

References:

- 1) Hara, Y. *et al.*: *DNA Res.*, 10, 129 (2003).
- 2) Koga, H. *et al.*: *DNA Res.*, 11, 293 (2004).

For research use only, not for diagnostic use

Manufactured by



ProteinExpress Co., Ltd. and



Kazusa DNA Research Institute.

Kazusa DNA Res.Inst.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619