



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-KIAA0338 Rat Homologue Gene Product (Type I Brain 4.1, Protein 4.1N) Polyclonal Antibody, Rabbit

Code	: PRX-PBR-1002
Format	: Protein A purified Rabbit IgG
Presentation	: Lyophilized, reconstitute with 0.2 ml of distilled water
Immunogen	: GST fusion-KIAA0338 Rat Homologue Gene Product (260 amino acids: L518-E777): LKEPNSKLIHRDRDWERERRLPSSPASPSPKGTPPEKASERA GLREGSEEKVKPPRPRAPESDTGDEDQDQERDAVFLKDNH LAIERKCSSITVSSTSSLEAEVDFTVIGDYHGGAFEDFSRSLP ELDRDKSDSETEGLVFARDLKGPSSQEDES GGIEDSPDRG ACSTPELPQFESVKAETMTVSSLAIRKKIEPEAMLQSRVSTA DSTQVDGGAPAAKDFMTTPPCITTETISTTMENSLKSGKGA AAMIPGPQTVATE
Buffer	: 0.15 M NaCl, 0.01 M Phosphate Buffer
Preservatives	: 0.05 % NaN ₃
Quantity	: 0.1 mg / vial
Applications	: ELISA 1:1000 Immunoblotting 1:1000 Immunohistochemistry 1:25
Storage	: 4°C

For research use only, not for diagnostic use

Manufactured by ProteinExpress Co., Ltd.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619