



COSMO BIO CO., LTD.
Inspiration for Life Science

Anti-KIAA0302 Rat Homologue Gene Product (β SpIII) Polyclonal Antibody, Rabbit

Code	: PRX-PBR-1003
Format	: Protein A purified Rabbit IgG
Presentation	: Lyophilized, reconstitute with 0.2 ml of distilled water
Immunogen	: GST fusion-KIAA0302 Rat Homologue Gene Product (β SpIII)(292 amino acids: Q2097-K2388): QPPTSEPMASQPEGSLVDGQRVLDTAWDGTQSKLPPSTQA PSINGVCTDTESSQPLLEQQRLEQSNVPEGPGSGTGDESS GPRGERQTLPRGPAPSPMPQSRSSSAHVATLPARGAELS AQEQMEGTLCRKQEMEAFNKKAAANRSWQNVYCVLRGSL GFYKDARAASAGVPYHGEVPVSLARAQGSVAFDYRK RKHV FKLGLQDGKEYLFQAKDEAEMSSWLRVVNAAIATASSASGE PEEPVVPSASRGLTRAMTPPVVSQPEGSIVLRSKDGRERE REKRFSFFKKNK
Buffer	: 0.01 M Phosphate Buffer containing 0.15 M NaCl, pH 7.4
Preservatives	: 0.05 % NaN_3
Quantity	: 0.1 mg / vial
Applications	: Immunohistochemistry 1:5
Storage	: 4°C

For research use only, not for diagnostic use

Manufactured by ProteinExpress Co., Ltd.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

URL: <http://www.cosmobio.co.jp>

e-mail: export@cosmobio.co.jp

[Outside Japan] Phone : +81-3-5632-9617

[国内連絡先] Phone : +81-3-5632-9610

FAX : +81-3-5632-9618

FAX : +81-3-5632-9619