



SPECIFICATION SHEET

**BIOTINYLATED ANTI-HUMAN β 2-MICROGLOBULIN
MOUSE MONOCLONAL ANTIBODY**

CATALOG NUMBER: IDL034

CLONE: EMR B6-12 (Biotin)

LOT NO:

QUANTITY: 50 μ g IgG / 100uL

SPECIFICITY: RECOGNIZE HUMAN β 2-MICROGLOBULIN

IMMUNOGEN: RECOMBINANT HUMAN β 2-MICROGLOBULIN

ISOTYPE: IgG1

FORMAT: PBS (pH7.0), 0.01%NaN₃

CONJUGATION: BIOTIN

APPLICATIONS: IMMUNOHISTOCHEMISTRY; 1:10-1:20
WESTERN BLOTING; 1:100-1:200

STORAGE: STORE AT 2~8°C.
STORE BELOW -20°C FOR LONG-TERM STORGE

REFERENCE: 1) KOMORI H, *et al.*, Clin. Cancer Res. 12(9), 2689-2697 (2006)

FOR RESEARCH USE ONLY, NOT FOR USE IN DIAGNOSTIC PROCEDURES

