



## **Anti Human Interferon Beta Monoclonal Antibody (Mouse)**

<b>Catalog No.:</b>	<b>SIM-2ZCKIB2</b>
<b>Clone No.:</b>	<b>Hyb-4020</b> <b>The splenic lymphocytes from BALB/c mouse immunized with E. coli expression ricombinant Interferon beta were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-4020) were screened.</b>
<b>Source:</b>	<b>BALB/C Mouse implantation ascites</b>
<b>Specificity:</b>	<b>This antibody specifically reacts with Human Interferon Beta</b>
<b>Purify:</b>	<b>Ammonium sulfate precipitation</b>
<b>Subclass:</b>	<b>IgG1</b>
<b>Buffer:</b>	<b>Physiological salt solution containing 0.1% Sodium Azide</b>
<b>Quantity:</b>	<b>0.5mL [0.5mg/mL]</b>
<b>Application:</b>	<b>ELISA</b>
<b>Storage:</b>	<b>Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.</b>

***For research use only; not for use as a diagnostic.***





コスモ・バイオ株式会社

SIM #2ZCKIB2

抗ヒト インターフェロン モノクローナル抗体(マウス)  
Anti IFN  $\beta$

コ ー ド : 2ZCKIB2

クローン番号 : Hyb- 4020

由 来 : 大腸菌発現リコンビナント - ヒトインターフェロン を免疫した  
BALB/c マウスの脾臓細胞とマウスミエローマ細胞 (NS-1) を  
融合して得たハイブリドーマを BALB/c マウスに接種して  
得られた腹水

特 異 性 : ヒトインターフェロン に特異的に反応

精 製 : 硫酸アンモニウム塩析精製

サブクラス : IgG1

溶 媒 : 0.1%アジ化ナトリウムを含む生理的食塩水溶液

蛋 白 濃 度 : 0.5 mg/mL

容 量 : 0.5 mL

貯 法 及 び : 2~8 3 ヶ月間

有 効 期 間 : -20 以下 3 年間

(凍結融解の繰り返しは避けてください)

販売元:



コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル  
TEL: (03) 5632-9610 FAX: (03) 5632-9619  
e-mail: mail@cosmobio.co.jp  
URL: <http://www.cosmobio.co.jp/>