



Anti Human Interferon Gamma Monoclonal Antibody (Mouse)

Catalog No.:	SIM-2ZCKIG2
Clone No.:	Hyb-6008 The splenic lymphocytes from BALB/c mouse immunized With E. coli expression ricombinant Human Interferon Gamma were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-6008) were screened.
Source:	BALB/C Mouse implantation ascites
Specificity:	This antibody specifically reacts with Human Interferon Gamma
Purify:	Ammonium sulfate precipitation
Subclass:	IgG1
Buffer:	Physiological salt solution containing 0.1% Sodium Azide
Quantity:	0.5mL [0.5mg/mL]
Application:	ELISA
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.





コスモ・バイオ株式会社

SIM #2ZCKIG2

抗ヒト インターフェロン モノクローナル抗体(マウス)
Anti IFN γ

コ ー ド : 2ZCKIG2

クローン番号 : Hyb- 6008

由 来 : 大腸菌発現リコンビナント - ヒトインターフェロン を
免疫した BALB/c マウスの脾臓細胞とマウスミエローム細胞
(NS-1) を融合して得たハイブリドーマを BALB/c マウスに
接種して得られた腹水

特 異 性 : ヒトインターフェロン に特異的に反応

精 製 : 硫酸アンモニウム塩析精製

サブクラス : IgG1

溶 媒 : 0.1%アジ化ナトリウムを含む生理的食塩水溶液

蛋 白 濃 度 : 0.5 mg/mL

容 量 : 0.5 mL

貯 法 及 び : 2~8 3 ヶ月間

有 効 期 間 : -20 以下 3 年間

(凍結融解の繰り返しは避けてください)

販売元:



コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
TEL: (03) 5632-9610 FAX: (03) 5632-9619
e-mail: mail@cosmobio.co.jp
URL: <http://www.cosmobio.co.jp/>