



Anti Human Blood group B agglutinin Monoclonal Antibody (Mouse)

Catalog No.:	SIM-2ZHBGB1
Clone No.:	Hyb-5362B The splenic lymphocytes from BALB/c mouse immunized With Human erythrocyte were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-5632B) were screened.
Source:	BALB/C Mouse implantation ascites
Specificity:	This antibody specifically reacts with Human blood group B agglutinin
Purify:	Ammonium sulfate precipitation
Subclass:	IgM
Buffer:	Physiological salt solution containing 0.1% Sodium Azide
Quantity:	0.5mL [0.5mg/mL]
Application:	Agglutinating
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.





コスモ・バイオ株式会社

SIM #2ZHBGB1

抗ヒト血液型 B モノクローナル抗体(マウス)
Anti Blood Group B

コード : 2ZHBGB1

クローン番号 : Hyb- 5362B

由来 : ヒト赤血球を免疫した BALB/c マウスの脾臓細胞とマウス
ミエローマ細胞 (NS-1) を融合して得たハイブリドーマを
BALB/c マウスに接種して得られた腹水

特異性 : ヒト B 型赤血球に対し凝集活性を示す

精製 : 硫酸アンモニウム塩析精製

サブクラス : IgM

溶媒 : 0.1%アジ化ナトリウムを含む生理的食塩水溶液

蛋白濃度 : 0.5 mg/mL

容量 : 0.5 mL

貯法及び : 2~8°C 3ヶ月間

有効期間 : -20°C以下 3年間

(凍結融解の繰り返しは避けてください)

販売元:



コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
TEL: (03) 5632-9610 FAX: (03) 5632-9619
e-mail: mail@cosmobio.co.jp
URL: <http://www.cosmobio.co.jp/>