



Anti Human Blood group Lea agglutinin Monoclonal Antibody (Mouse)

Catalog No.:	SIM-2ZHBGLA
Clone No.:	Hyb-4620 The splenic lymphocytes from BALB/c mouse immunized With Human saliva origin blood group substance were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-4620) were screened.
Source:	BALB/C Mouse implantation ascites
Specificity:	This antibody specifically reacts with Human blood group Lea agglutinin
Purify:	Ammonium sulfate precipitation
Subclass:	IgG1
Buffer:	Physiological salt solution containing 0.1% Sodium Azide
Quantity:	0.5mL [0.5mg/mL]
Application:	Agglutinating
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.





コスモ・バイオ株式会社

SIM #2ZHBGLA

抗ヒト血液型 Lea モノクローナル抗体(マウス)
Anti Blood Group Lea

コード : 2ZHBGLA

クローン番号 : Hyb- 4620

由来 : ヒト唾液由来血液型物質を免疫した BALB/c マウスの脾臓細胞とマウスミエローマ細胞 (NS-1) を融合して得たハイブリドーマを BALB/c マウスに接種して得られた腹水

特異性 : ヒト Le(a+b-) 型赤血球に対し凝集活性を示す (クームス血清存在下)

精製 : 硫酸アンモニウム塩析精製

サブクラス : IgG1

溶媒 : 0.1%アジ化ナトリウムを含む生理的食塩水溶液

蛋白濃度 : 0.5 mg/mL

容量 : 0.5 mL

貯法及び : 2~8°C 3ヶ月間

有効期間 : -20°C以下 3年間

(凍結融解の繰り返しは避けてください)

販売元:



コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
TEL: (03) 5632-9610 FAX: (03) 5632-9619
e-mail: mail@cosmobio.co.jp
URL: <http://www.cosmobio.co.jp/>