



## **Anti Human Glycophorin B Monoclonal Antibody (Mouse)**

<b>Catalog No.:</b>	<b>SIM-2ZHBGPB</b>
<b>Clone No.:</b>	<b>Hyb-5362C</b> <b>The splenic lymphocytes from BALB/c mouse immunized With Human erythrocyte were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-5632C) were screened.</b>
<b>Source:</b>	<b>BALB/C Mouse implantation ascites</b>
<b>Specificity:</b>	<b>This antibody specifically reacts with Human Glycophorin B (This antibody coheres with all of ABO blood group and MN blood group)</b>
<b>Purify:</b>	<b>Ammonium sulfate precipitation</b>
<b>Subclass:</b>	<b>IgG3</b>
<b>Buffer:</b>	<b>Physiological salt solution containing 0.1% Sodium Azide</b>
<b>Quantity:</b>	<b>0.5mL [0.5mg/mL]</b>
<b>Application:</b>	<b>ELISA</b>
<b>Storage:</b>	<b>Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles.</b>

***For research use only; not for use as a diagnostic.***





コスモ・バイオ株式会社

SIM #2ZHBGPB

抗ヒト グリコホリン B モノクローナル抗体(マウス)  
Anti Glycophorin B

コード : 2ZHBGPB

クローン番号 : Hyb- 5362C

由来 : ヒト赤血球を免疫した BALB/c マウスの脾臓細胞とマウス  
ミエローマ細胞 (NS-1) を融合して得たハイブリドーマを  
BALB/c マウスに接種して得られた腹水

特異性 : ヒトグリコホリン B に特異的に反応  
( ABO 型、MN 型全てと凝集 )

精製 : 硫酸アンモニウム塩析精製

サブクラス : IgG3

溶媒 : 0.1%アジ化ナトリウムを含む生理的食塩水溶液

蛋白濃度 : 0.5 mg/mL

容量 : 0.5 mL

貯法及び : 2~8 3ヶ月間

有効期間 : -20 以下 3年間

(凍結融解の繰り返しは避けてください)

販売元:



コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル  
TEL: (03) 5632-9610 FAX: (03) 5632-9619  
e-mail: mail@cosmobio.co.jp  
URL: <http://www.cosmobio.co.jp/>