



Anti HBc Monoclonal Antibody Antigenic determinant alpha (Mouse)

Catalog No.:	SIM-2ZHC21
Clone No.:	Hyb-3105 The splenic lymphocytes from BALB/c mouse immunized With purified HBc antigen were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-3105) were screened.
Source:	BALB/C Mouse implantation ascites
Specificity:	This antibody specifically reacts with HBc antigen (Antigenic determinant alpha)
Purify:	Ammonium sulfate precipitation
Subclass:	IgG3
Buffer:	Physiological salt solution containing 0.1% Sodium Azide
Quantity:	0.5mL [0.5mg/mL]
Application:	Western Blot, ELISA
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.





コスモ・バイオ株式会社

SIM #2ZHC21

抗 HBc モノクローナル抗体 抗原決定基 (マウス)
Anti HBc determinant α

コード : 2ZHC21

クローン番号 : Hyb- 3105

由来 : 精製 HBc 抗原を免疫した BALB/c マウスの脾臓細胞とマウス
ミエローム細胞 (NS-1) を融合して得たハイブリドーマを
BALB/c マウスに接種して得られた腹水

特異性 : HBc 抗原 (抗原決定基) に特異的に反応

精製 : 硫酸アンモニウム塩析精製

サブクラス : IgG3

溶媒 : 0.1%アジ化ナトリウムを含む生理的食塩水溶液

蛋白濃度 : 0.5 mg/mL

容量 : 0.5 mL

貯法及び : 2~8 3ヶ月間

有効期間 : -20 以下 3年間

(凍結融解の繰り返しは避けてください)

販売元:



コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
TEL: (03) 5632-9610 FAX: (03) 5632-9619
e-mail: mail@cosmobio.co.jp
URL: http://www.cosmobio.co.jp/