



Anti Human secretor IgA Monoclonal Antibody (Mouse)

Catalog No.:	SIM-2ZHIGSA
Clone No.:	Hyb-5831 The splenic lymphocytes from BALB/c mouse immunized with human colostrum origin secretor IgA were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-5831) were screened.
Source:	BALB/C Mouse implantation ascites
Specificity:	This antibody specifically reacts with Human secretor IgA
Purify:	Ammonium sulfate precipitation
Subclass:	IgG1
Buffer:	Physiological salt solution containing 0.1% Sodium Azide
Quantity:	0.5mL [0.5mg/mL]
Application:	ELISA
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.





コスモ・バイオ株式会社

SIM #2ZHIGSA

抗ヒト分泌型 IgA モノクローナル抗体(マウス)
Anti IgA Secretory

コード : 2ZHIGSA

クローン番号 : Hyb- 5831

由来 : ヒト初乳由来分泌型 IgA を免疫した BALB/c マウスの脾臓細胞と
マウスミエローマ細胞 (NS-1) を融合して得たハイブリドーマを
BALB/c マウスに接種して得られた腹水

特異性 : ヒト分泌型 IgA に特異的に反応

精製 : 硫酸アンモニウム塩析精製

サブクラス : IgG1

溶媒 : 0.1%アジ化ナトリウムを含む生理的食塩水溶液

蛋白濃度 : 0.5 mg/mL

容量 : 0.5 mL

貯法及び : 2~8 3ヶ月間

有効期間 : -20 以下 3年間

(凍結融解の繰り返しは避けてください)

販売元:



コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
TEL : (03) 5632-9610 FAX : (03) 5632-9619
e-mail : mail@cosmobio.co.jp
URL : <http://www.cosmobio.co.jp/>