



Anti HP Urease B Monoclonal Antibody (Mouse)

Catalog No.:	SIM-2ZHPUB
Clone No.:	Hyb-K2031B The splenic lymphocytes from BALB/c mouse immunized with <i>Helicobacter Pylori</i> (HP) origin Urease were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-K2031B) were screened.
Source:	BALB/C Mouse implantation ascites
Specificity:	This antibody specifically reacts with HP Urease large subunit (UreB)
Purify:	Ammonium sulfate precipitation
Subclass:	IgG1
Buffer:	Physiological salt solution containing 0.1% Sodium Azide
Quantity:	0.5mL [0.5mg/mL]
Application:	Western Blot, IHC frozen section, ELISA
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.





コスモ・バイオ株式会社

SIM #2ZHPUB

抗 HP ウレアーゼ B モノクローナル抗体(マウス)
Anti HP Urease B

コード : 2ZHPUB

クローン番号 : Hyb- K2031B

由来 : Helicobacter Pylori (HP) 菌体由来ウレアーゼを免疫した BALB/c
マウスの脾臓細胞とマウスミエローマ細胞 (NS-1) を融合して
得たハイブリドーマを BALB/c マウスに接種して得られた腹水

特異性 : HP ウレアーゼ large subunit (UreB) に特異的に反応

精製 : 硫酸アンモニウム塩析精製

サブクラス : IgG1

溶媒 : 0.1%アジ化ナトリウムを含む生理的食塩水溶液

蛋白濃度 : 0.5 mg/mL

容量 : 0.5 mL

貯法及び : 2~8 3ヶ月間

有効期間 : -20 以下 3年間

(凍結融解の繰り返しは避けてください)

販売元:



コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
TEL: (03) 5632-9610 FAX: (03) 5632-9619
e-mail: mail@cosmobio.co.jp
URL: http://www.cosmobio.co.jp/