



Anti Mouse RBP-Jk Monoclonal Antibody (Rat)

Catalog No.:	SIM-2ZRBP1
Clone No.:	Hyb-K0043 The splenic lymphocytes from PVG/cSea rat immunized with baculovirus expression Mouse RBP-Jk were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-K0043) were screened.
Source:	BALB/C Mouse implantation ascites
Specificity:	This antibody specifically reacts with RBP-Jk
Purify:	Ammonium sulfate precipitation
Subclass:	IgG1
Buffer:	Physiological salt solution containing 0.1% Sodium Azide
Quantity:	0.5mL [0.5mg/mL]
Application:	Immunohistochemistry, Immunoprecipitation, ELISA
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

