



Anti Mouse RBP-Jk Monoclonal Antibody (Rat)

Catalog No.:	SIM-2ZRBP2
Clone No.:	Hyb-T6709 The splenic lymphocytes from PVG/cSea rat immunized with E.coli expression Mouse RBP-Jk were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-T6709) were screened.
Source:	Culture supernatant (culture medium : ASF-104)
Specificity:	This antibody specifically reacts with RBP-Jk
Purify:	Centrifugation
Subclass:	IgG2a
Buffer:	Culture supernatant containing 0.1% Sodium Azide
Quantity:	2mL
Application:	Western Blot, Immunohistochemistry, ELISA
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

