



Anti Mouse RBP-Jk Monoclonal Antibody (Rat)

Catalog No.:	SIM-2ZRBP2
Clone No.:	Hyb-T6709 The splenic lymphocytes from PVG/cSea rat immunized with E.coli expression Mouse RBP-Jk were fused to mouse myeloma NS-1 cells and the hybrid cells (Hyb-T6709) were screened.
Source:	Culture supernatant (culture medium : ASF-104)
Specificity:	This antibody specifically reacts with RBP-Jk
Purify:	Centrifugation
Subclass:	IgG2a
Buffer:	Culture supernatant containing 0.1% Sodium Azide
Quantity:	2mL
Application:	Western Blot, Immunohistochemistry, ELISA
Storage:	Maintain at 4-8 degrees (-20 degrees for long term storage) in undiluted aliquots.Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.





コスモ・バイオ株式会社

SIM #2ZRBP2

抗マウス RBP-J κ モノクローナル抗体(ラット)
Anti RBP-J κ

コード : 2ZRBP2

クローン番号 : Hyb- T6709

由来 : 大腸菌発現マウス RBP-J κ を免疫した PVG/cSea ラットの脾臓細胞とマウスミエローマ細胞 (NS-1) を融合して得たハイブリドーマを培養して得られた培養上清

特異性 : RBP-J κ に特異的に反応

精製 : ゴミ取り遠心精製

サブクラス : IgG2a

溶媒 : 0.1%アジ化ナトリウムを含む培養上清 (培養液 : ASF-104)

蛋白濃度 :

容量 : 2 mL

貯法及び : 2~8°C 3ヶ月間

有効期間 : -20°C以下 3年間

(凍結融解の繰り返しは避けてください)

販売元:



コスモ・バイオ株式会社

〒135-0016 東京都江東区東陽 2-2-20 東陽駅前ビル
TEL : (03) 5632-9610 FAX : (03) 5632-9619
e-mail : mail@cosmobio.co.jp
URL : <http://www.cosmobio.co.jp/>