



Antibody to Human B-type Natriuretic Peptide

Abbreviation : RP BNP

Code No. : 7624

Price : 50 μ l ¥50,000

Quality : Rabbit serum

Solution : Serum, Saline (50 μ l)

Character :		Reactivity
to	Human BNP (1-32)	100 %
	Human ANP (1-28)	< 0.04 %
	Human CNP (1-22)	< 0.04 %

Example :	Add successively.	*Assay buffer
	100 μ l of sample in the assay buffer*	100 mM phosphate buffer pH 7.4
	100 μ l of $\times 1000$ diluted antiserum	1 mM EDTA dipotassium
	Incubate at 4°C for 24 hours.	10 mM glycine
	Add 100 μ l of 125 I-HBNP solution.	10 mM ϵ -aminocaproic acid
	Incubate at 4°C for 24 hours.	1 mM sodium azide
	Separate by the double antibody method.	0.6 % triton X-100

Storage & Handling : The solution should be stored below -20°C.

Repeated freezing and thawing should be avoided.

Note : Research use only.

Reference : Togashi, K. et al. Jpn. J. Clin. Pathol. **39**, 283-288 (1991)

