




Monoclonal Antibody to Human PTH C-terminal

Code No.: 7648
Quantity: 200 µl (200 µg)
Price: ¥40,000
Clone: 5A7
Quality: BALB/c mouse IgG₁ κ (Purified by protein A chromatography)
Solution: 1mg/ml IgG, 0.1%NaN₃, PBS, pH7.4
Character: This antibody is specific to C-terminal (69-84) of human PTH.



Storage: The solution should be stored below -20°C.
Repeated freezing and thawing should be avoided.

Note: Research use only.