




Monoclonal Antibody to Rat PTH N-terminal

Code No.: 7667
Quantity: 200 μ l(200 μ g)
Price: ¥40,000
Clone: 8A2
Quality: BALB/c mouse IgG₁ κ (Purified by protein A chromatography)
Solution: 1mg/ml IgG, 0.1%NaN₃, PBS, pH7.4
Character: This antibody is specific to N-terminal (2-34) of rat PTH.



Storage: The solution should be stored below -20°C.
Repeated freezing and thawing should be avoided.

Note : Research use only.