



Monoclonal Mouse IgE anti-DNP

Abbreviation : M-IgE

Code No. : 7676

Price : 100 μ l (100 μ g) ¥35,000

Clone : 2-9(4)

Quality : Purified from ascites by column chromatographies.

Solution : 1 mg/ml IgE, 0.1 % NaN₃, PBS, pH 7.4 (100 μ l)

Character : This preparation shows single band on SDS-PAGE.
This antibody is specific for DNP.

Storage & Handling : The solution should be stored below -20°C.

Handling : Repeated freezing and thawing should be avoided.

Note : Research use only.

References : 1) Saito, H. et al. J. Immunol. **138**, 3927 (1987)
2) Nakamura, H. et.al. Jpn. J. Allergol. **38**(12), 1359 (1989)