




Monoclonal Antibody to Human Midkine

Code No. : 80037
Product name: MAB TO HUMAN MIDKINE (IP-10)
Quantity: 200 μ L (200 μ g)
Clone: IP-10
Quality: Mouse IgG1 (Purified by protein A chromatography)
Solution: 1 mg/mL IgG, 0.1% NaN_3 , PBS, pH 7.4 (200 μ L)
Character:



Storage: The solution should be stored below -20°C .
Repeated freezing and thawing should be avoided.
Note: Research use only.
Reference: Tokuyama W, et al., Pathol Int. 57, 260-267(2007).