




Monoclonal Antibody to Human Midkine

Code No. : 80038
Product name: MAB TO HUMAN MIDKINE (SC-2)
Quantity : 200 μ L (200 μ g)
Clone : SC-2
Quality : Mouse IgG1 (Purified by protein A chromatography)
Solution : 1 mg/mL IgG, 0.1% NaN_3 , PBS, pH 7.4 (200 μ L)
Character:



Storage: The solution should be stored below -20°C .
Repeated freezing and thawing should be avoided.

Note: Research use only.