




Monoclonal Antibody to Human Midkine

Code No. : 80039
Product name: MAB TO HUMAN MIDKINE (SC-4)
Quantity: 200 μ L (200 μ g)
Clone: SC-4
Quality: Mouse IgG1 (Purified by protein A chromatography)
Solution: 1 mg/mL IgG, 0.1% NaN₃, PBS, pH 7.4 (200 μ L)
Character:



Storage: The solution should be stored below -20°C.
Repeated freezing and thawing should be avoided.
Note: Research use only.