




Monoclonal Antibody to Lewis^X

Code No. : 80040
Product name: MAB TO LEWIS X
Quantity: 200 μ L
Clone: 4C9
Quality: Rat IgM
Solution: Ascites (200 μ L)
Character: Working dilution is about 1: 100



Storage: The solution should be stored below -20°C .
Repeated freezing and thawing should be avoided.

Note: Research use only.

Reference: 1) Yoshinaga, K. et al., Differentiation. 48, 75-82 (1991)
2) Matsumoto, H. et al., Cancer. 69, 2084-90 (1992)