



FITC Conjugated Monoclonal Antibody to Nephritis-associated plasmin receptor

Code No.: 80041

Product name: MAB TO NEPHRITIS-ASSOCIATED PLASMIN-R

Quantity: 200 μ L (200 μ g protein)

Clone: 1F10

Quality: FITC conjugated BALB/c IgG2 κ

Solution: 1mg/mL protein in PBS, pH7.6

Storage: The solution should be stored below -20°C.
Repeated freezing and thawing should be avoided.

Note: Research use only.

Reference: Oda T, et al., J Am Soc Nephrol. 2005 Jan; 16(1):247-54.
Yoshizawa N, et al., J Am Soc Nephrol. 2004 Jul; 15(7):1785-93.
Yamakami K, et al., Methods. 2000 Jun; 21(2):185-97.