

Human Growth Hormone

Purified Factor Product Number HC88851 100 ug

DESCRIPTION:

Human Growth Hormone, fully biologically active.

SOURCE: E. Coli

Purity > 95% by HPLC analysis

BIOLOGICAL ACTIVITY:

Fully biologically active when compared to standards.

Measured at 2.6 units per mg based upon WHO reference standard.

HANDLING AND STORAGE:

Purified HGH factor is provided as semi-frozen liquid, in PBS. Store FROZEN at -20° C. dilute in PBS or other buffered solutions. Store small aliquots at -20° C for long-term storage, at high concentration.

Avoid repeat freeze-thaw.

Before opening vial to reconstitute, it is recommended to microcentrifuge vial for 30 seconds -

WARRANTY:

Products sold hereunder are warranted only to conform to the quantity and contents stated on the label at the time of delivery to the customer. There are no warranties, expressed or implied, which extend beyond the description on the label of the product.

RESEARCH USE ONLY -NOT for DRUG USE