

AAA

#RHF290

**ANTIGENIX AMERICA INC.**

**1-800-558-1008 Phone**

**1-631-424-0411 Fax**

**PRODUCT INFORMATION SUMMARY**

**Anti-EGF**

**(Polyclonal)**

Purified Polyclonal Antibody Product Number RHF290 0.1 mg

Biotin Conjugate Product Number RHF290-B 50 ug

**\*Antigen-Affinity Purified**

**SPECIFICITY:**

Rabbit polyclonal antibodies to human Epidermal growth Factor(EGF), produced from highly purified recombinant human EGF immunogen. Antiserum fraction further purified by **Antigen-Affinity** chromatography to achieve superior titre and specificity.

**SOURCE:** Rabbit antiserum.

**RESEARCH APPLICATIONS:**

Identification of human EGF in body fluids and tissue sections.

**ELISA:** Antibody concentration of 0.5 ug/mL (100 uL/well) was reacted in ELISA to detect 0.4 ng/well recombinant EGF.

**Biotin conjugate:** use @ 0.2 - 0.3 ug/mL

**NEUTRALIZING ACTIVITY:** Antibody concentration of 0.35 ug/mL produced one-

half maximal inhibition of 2.0 ng/mL human EGF.

**WESTERN BLOT:** Antibody concentration of 0.2 ug/mL was used to detect 1.5 ng/lane recombinant human EGF under both reducing and non-reducing conditions.

**HANDLING AND STORAGE:**

Purified polyclonal antibody provided lyophilized from PBS, preservative-free. Reconstitute with 100 uL sterile or distilled water. **For Biotin Conjugate:** Reconstitute with 100 uL sterile water **containing 0.1% BSA.** These preparations should be diluted in a protein-containing or other stabilizing medium to a concentration suitable for use in specific protocols. Small aliquots should be frozen at -20° C for long-term storage. Reagents in a liquid state are stable for approximately one month when stored at 2-8° C. Avoid repeat freeze-thaw cycles.

**WARRANTY:**

Products sold hereunder are warranted only to conform to the quantity and contents stated on the label at the time of delivery to the customer. There are no warranties, expressed or implied, which extend beyond the product label description.

**RESEARCH USE ONLY**

**NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE.**