



AAA \*RMF600

## PRODUCT INFORMATION SUMMARY

### Anti-Mouse IL-6 (Polyclonal)

Purified Polyclonal Antibody      Product Number      RMF600      0.5 mg

#### **SPECIFICITY:**

Rabbit polyclonal antibodies to mouse Interleukin 6 (IL-6), produced from highly purified recombinant mouse IL-6 immunogen. Antiserum fraction purified by protein A chromatography.

**SOURCE:**              Rabbit antiserum.

#### **RESEARCH APPLICATIONS:**

Identification of mouse IL-6 in body fluids and tissue sections.

ELISA: 1:8,000 dilution from 1mg/mL antibody concentration was reacted in direct ELISA to detect 1 ug/mL murine IL-6 coated on microplate.

**NEUTRALIZING ACTIVITY:** Not tested

#### **HANDLING AND STORAGE:**

Purified polyclonal antibody provided lyophilized from PBS, preservative-free. Reconstitute with 500 uL distilled water. These preparations should be diluted in buffered solutions (PBS or similar) or other stabilizing medium to a concentration suitable for use in specific protocols. Small aliquots should be frozen at -20° C for long-term storage, or preservative added as suitable for application. All reagents in a liquid state should be stored at 2-8° C when not in use, with addition of preservative (0.1% sodium azide or similar) to store in liquid state.. Avoid repeat freeze-thaw cycles.

#### **WARRANTY:**

Products sold hereunder are warranted only to conform to the quantity and contents stated on the label at the time of delivery to the customer. There are no warranties, expressed or implied, which extend beyond the product label description.

**RESEARCH USE ONLY  
NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE.**