



Product Data Sheet for:  
**Sheep anti-human Neutrophil Elastase**  
**Affinity purified IgG**

**Catalogue #:** SANE-AP

**Lot #:** XXXX

**Expiry:** XXXX

Rev: 103

**Immunogen:** Neutrophil elastase (active-site blocked) prepared from human sputum.

**Presentation:** 0.127 ml of immunoglobulin affinity purified from sheep antiserum on immobilized elastase (active-site blocked). The IgG concentration is 3.94 mg/ml in 10 mM HEPES, pH 7.4, containing 0.15 M NaCl and 50% glycerol (v/v).

**Storage:** Store between -10 and -20°C. Gently warm and mix before use.

**Specificity:** This antibody is specific for elastase as demonstrated by immunoelectrophoresis and ELISA.

Related products:

Product: <b>SANE-IG</b>	Sheep anti-human Neutrophil Elastase, purified whole IgG from antiserum	10 mg vial
Product: <b>SANE-HRP</b>	Sheep anti-human Neutrophil Elastase, peroxidase labelled whole IgG	0.2 mg vial

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

**For Research Use Only**