



Product Data Sheet for:
**Sheep anti-human Neutrophil Elastase
Purified IgG**

Catalogue #: SANE-IG
Lot #: XXXX
Expiry: XXXX

Rev: 1003

Immunogen: Neutrophil elastase (active-site blocked) prepared from human sputum.

Presentation: 1 ml of immunoglobulin purified from sheep antiserum. The IgG concentration is 10 mg/ml in 10 mM HEPES, pH 7.4, containing 0.15 M NaCl and 50% glycerol (v/v).

Storage: Store in the freezer between -10 and -20°C. Gently mix before use.

Specificity: IgG is monospecific for elastase when analyzed by immunoelectrophoresis and ELISA.

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

For Research Use Only