



Product Data Sheet for:
Sheep anti-Human Vimentin
Affinity Purified IgG

Catalogue #: SAVM-AP

Lot #: XXXX

Expiry: XXXX

Rev: 103

Immunogen: Purified recombinant vimentin (Vm) produced in bacteria. This protein constitutes a 20 kDa fragment from the amino-terminus of vimentin.

Presentation: 0.25 ml of immunoglobulin from sheep antiserum affinity purified on recombinant vimentin fragment immobilized on agarose beads. The IgG concentration is 2.0 mg/ml in 10mM HEPES, pH 7.4, 0.15M NaCl and 50% glycerol (v/v). This product contains no preservative.

Storage: Store in the freezer between -10 to -20°C. Gently mix before use.

Specificity: This antibody is specific for vimentin as demonstrated by immunoblotting and immunohistochemistry. No reactivity was seen with citrated platelet-poor plasma by immunoelectrophoresis.

Application: This product can be used for immunostaining and immunoblotting applications for human vimentin.

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory

For Research Use Only