


Affinity Biologicals™ INC.

 1395 Sandhill Drive. Ancaster, Ontario, Canada L9G 4V5
 905-304-9896 • 800-903-6020 • fax 905-304-9897

**** REPRESENTATIVE DATASHEET ****

 Product Data Sheet for:
**Sheep anti-Human Vimentin
 Purified IgG**

Catalogue #: SAVM-IG
Lot #: XXXXXXXX
Expiry: MM/YYYY

Rev: 504

Immunogen: Purified recombinant human vimentin (Vm) produced in bacteria. This protein constitutes a 20 kDa fragment from the amino-terminus of vimentin.

Presentation: XX ml of immunoglobulin purified from sheep antiserum. The IgG concentration is XX mg/ml in 10 mM HEPES, pH 7.4, 0.15 M NaCl and 50% glycerol (v/v). Total antibody is 10 mg.

Storage: Store in the freezer between -10 to -20°C. Gently mix before use.

Specificity: This antibody is specific for vimentin as demonstrated by immunoelectrophoresis and ELISA.

Application: This product can be used as a source of antibodies to human vimentin.

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory.

For Research Use Only
Related Products:

| | | |
|----------------------------|--|-------------|
| Product: SAVM-AP | Sheep anti-human Vimentin, affinity purified IgG (APIgG) | 0.5 mg vial |
| Product: SAVM-APHRP | Sheep anti-human Vimentin, peroxidase-conjugated APIgG | 0.1 mg vial |