



Affinity Biologicals™ INC.

1395 Sandhill Drive. Ancaster, Ontario, Canada L9G 4V5
905-304-9896 • 800-903-6020 • fax 905-304-9897

**** REPRESENTATIVE DATASHEET ****

Product Data Sheet for:

**Sheep anti-Human Vitronectin
Purified IgG**

Catalogue #: SAVN-IG
Lot #: XXXXXXXXX
Expiry: MM/YYYY

Rev: 302

Immunogen: Purified vitronectin (Vn) from human plasma. This protein is also known as S-Protein.

Presentation: XX ml of immunoglobulin from sheep antiserum. The IgG concentration is XX mg/ml in 0.01 M HEPES, pH 7.4, containing 0.1 M NaCl and 50% glycerol (v/v). Total antibody is 10 mg.

Storage: Store in the freezer between -10 to -20°C. Gently mix before use.

Specificity: This antibody is specific for vitronectin as demonstrated by immunoelectrophoresis and ELISA.

Application: This product can be used as a source of antibodies to human vitronectin.

Optimum conditions for use of this antibody in analytical procedures should be established in each laboratory.

For Research Use Only

Related products:

Product: SAVN-AP	Sheep anti-human Vitronectin, affinity purified IgG (APIgG)	0.5 mg vial
Product: SAVN-APHRP	Sheep anti-human Vitronectin, peroxidase-conjugated APIgG	0.1 mg vial