

## CryoSolutions 100%

*Akron Biotech Catalog# AK8899*

---

### Product Description

**CryoSolutions 100%** is a cryopreservation medium, sterile, filtered DMSO (Dimethyl Sulfoxide, USP/EP/JP) - highest purity grade >99.9%.

### Product Specifications:

- ✓ Contains >99.9% dimethyl sulfoxide (DMSO)
- ✓ Manufactured in compliance with cGMP guidelines
- ✓ DMSO meets Pharmacopeia standards (USP, EP, JP)
- ✓ Free of ions and heavy metals
- ✓ Tested for sterility and endotoxin

### Packaging:

- ✓ Packaged in vials made from borosilicate glass USP Type I
- ✓ Septum stoppers, crimp seal and flip-off caps – certified, sterilized—designed for aseptic operations.
- ✓ Sizes: 5 ml, 10 ml, 30 ml, 100 ml; Bulk sizes available

### Research Application:

**CryoSolutions 100%** is intended for cryopreservation of biological cells in cell culture, stem cells, mammalian cells and biological tissues.

**Storage Requirements:** Store at room temperature. Protect from light.

### Related Applications:

- ✓ Cell Culture
- ✓ Reagents and Supplements
- ✓ Cell Counting, Viability, and Cryopreservation
- ✓ Mammalian Cell Cultures
- ✓ Primary Cell Cultures

### Related Products:

CryoSolutions 10% (Cat. #AK8900)  
CryoSolutions MC (Cat. #AK8995)  
CryoSolutions DX (Cat. #AK8996)

CryoSolutions DX5 (Cat. #AK9544)  
AK-Freez (Cat. #AK8361)  
CryoSolve™ (Cat. #AK8217)