

AALTO



Purified Recombinant Streptolysin O (SLO) Antigen

PRODUCT CODE : AW 6066-b

SOURCE : Transformed E. Coli

STORAGE BUFFER : Tris-HCl buffered saline, pH 8.5

PURIFICATION : Metal affinity chromatography.

PROTEIN CONCENTRATION : 2 - 4 mg/ml (OD_{280nm} x 0.79)

TOTAL ANTIGENIC ACTIVITY : > 4000 IU/ml by internal indirect method

SPECIFIC ANTIGENIC ACTIVITY : > 1500 IU/mg protein.

PRESERVATIVE : None

STORAGE : Frozen (-70°C)
Avoid repeated freezing and thawing

Aalto Bio Reagents Ltd.

Church Lane, Rathfarnham Village
Dublin 14, Ireland.

Tel: +353-1-490 0685 Fax: +353-1-490 0122

E mail: info@aaltobioreagents.ie Web: www.aaltobioreagents.ie



1.8. EN ISO 9001:2000

