

AALTO



Purified and Stabilised Recombinant Streptolysin O (SLO) Antigen

PRODUCT CODE : AW 6093-b

SOURCE : Transformed E. Coli

STORAGE BUFFER : Tris-HCl buffered saline, pH 8.2

PURIFICATION : Metal affinity chromatography.

PROTEIN CONCENTRATION : 1-2 mg/ml (OD_{280nm} x 0.69)

TOTAL ANTIGEN ACTIVITY : >2500 IU/ml by latex enhanced immunoturbidimetric method

SPECIFIC ANTIGEN ACTIVITY : >1250 IU/mg protein.

PRESERVATIVE : < 0.1% sodium azide

STORAGE : Frozen (-70°C)
Avoid repeated freezing and thawing

Aalto Bio Reagents Ltd.

Church Lane, Rathfarnham Village
Dublin 14, Ireland.

Tel: +353-1-490 0685 Fax: +353-1-490 0122

E mail: info@aaltobioreagents.ie Web: www.aaltobioreagents.ie



1.8. EN ISO 9001:2000

