

# Anti-PFKFB3 (6-phosphofructo-2-kinase/fructose-2,6-biphosphatase 3) Polyclonal Antibody

Cat. #: 60B286

## Description:

The PFKFB3(6-phosphofructo-2-kinase/fructose-2,6-biphosphatase 3) is involved in synthesis and hydrolysis of fructose-2,6-bisphosphate. In the C-terminal section, PFKFB3 belongs to the phosphoglycerate mutase family.

## Immunogen/Specificity:

Polyclonal antibody produced in rabbits immunizing with a synthetic peptide corresponding to C-terminal residues of human PFKFB3(6-phosphofructo-2-kinase/fructose-2,6-biphosphatase 3)

## References

- Sakai,A., et al, J. Biochem. 119 (3), 506-511 (1996)  
Manzano,A., et al, Cytogenet. Cell Genet. 83 (3-4), 214-217 (1998)  
Chesney,J., et al, Proc. Natl. Acad. Sci. U.S.A. 96 (6), 3047-3052 (1999)  
Scanlan,M.J., et al, Int. J. Cancer 83 (4), 456-464 (1999)

Clone Number:

Isotype:

Species: Human

Storage and Stability: at -20oC

Storage buffer:

This antibody is stored in PBS, 0.01% sodium azide and 50% glycerol.

Preparation:

Purified by antigen-specific affinity chromatography.

Applications :

ELISA

Western Blotting (1µg/ml for 2hrs)