

## Anti-BRF1(Butyrate response factor 1) Polyclonal Antibody

Cat. #: 60B603

### Description:

BRF1(Butyrate response factor 1) is a probable regulatory protein involved in regulating the response to growth factors. BRF1 contains 2 C3H1-type zinc fingers.

### Immunogen/Specificity:

Polyclonal antibody produced in rabbits immunizing with a synthetic peptide corresponding to N-terminal residues of human BRF1(Butyrate response factor 1)

### References

Barnard,R.C.,et al, Nucleic Acids Res. 21 (15), 3580 (1993)  
Bustin,S.A., et al, DNA Cell Biol. 13 (5), 449-459 (1994)  
Ning,Z.Q.,et al, Eur. J. Immunol. 26 (10), 2356-2363 (1996)

Clone Number:

Isotype:

Species: Human

Storage and Stability: at -20oC

### Storage buffer:

This antibody is stored in PBS, 0.01% sodium azide and 50% glycerol.

### Preparation:

Purified by antigen-specific affinity chromatography.

### Applications :

ELISA

Western Blotting (1µg/ml for 2hrs)