

Anti-RAB8A(Ras-related protein Rab-8A)Polyclonal Antibody

Cat. #: 60B636

Description:

The RAB8A (Ras-related protein Rab-8A) may be involved in vesicular trafficking and neurotransmitter release. It interacts with MAP4K2 and SYTL4. The protein belongs to the small GTPase superfamily. Rab family.

Immunogen/Specificity:

Polyclonal antibody produced in rabbits immunizing with a synthetic peptide corresponding to N-terminal residues of human RAB8A (Ras-related protein Rab-8A)

References

Nimmo, E.R., et al, *Oncogene* 6 (8), 1347-1351 (1991)
Joberty, G., et al, *FEBS Lett.* 330 (3), 323-328 (1993)

Clone Number:

Isotype:

Species: Human

Storage and Stability: at -20°C

Storage buffer:

This antibody is stored in PBS, 0.01% sodium azide and 50% glycerol.

Preparation:

Purified by antigen-specific affinity chromatography.

Applications :

ELISA

Western Blotting (1 µg/ml for 2hrs)