

PTEN (Ab-380/382/383) Antibody

Cat. #: ANTY021056

Species: Human, Mouse and Rat

Quantity: 100ug

Concentration: 100ug/100ul

Storage and Stability:

Store at -20°C/1 year

Immunogen:

The antiserum was produced against synthesized non-phosphopeptide derived from human PTEN around the phosphorylation site of serine 380 and threonine 382/383 (R-Y-SP-D-TP-TP-D-S).

Specificity:

PTEN (Ab-380/382/383) antibody detects endogenous levels of total PTEN protein.

Tested Application: WB and IHC

Application Notes: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Raised In: Rabbit


Purity: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Storage buffer: WB: 1:500~1:1000 IHC:1:50~100

Form: Liquid

Reference:

Al-Khouri AM, et al. (2005) J Biol Chem. 280(42):35195-35202. Torres J, et al. (2001) J Biol Chem. 276(2): 993-998. Vazquez F, et al. (2000) Mol Cell Biol. 20(14): 5010-5018.



For research use only

Antagene, Inc. Toll Free: 1-866-964-2589, Fax: 1-888-225-1868, info@antageneinc.com,
www.antageneinc.com