

HDAC4 (Ab-632) Antibody

Cat. #: ANTY021141

Species: Human, Mouse and Rat

Quantity: 100ug

Concentration: 100ug/100ul

Storage and Stability:

Store at -20°C/1 year

Immunogen:

The antiserum was produced against synthesized non-phosphopeptide derived from human HDAC4 around the phosphorylation site of serine 632 (A-Q-SP-S-P).

Specificity:

HDAC4 (Ab-632) antibody detects endogenous levels of total HDAC4 protein.

Tested Application: WB

Application Notes: Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Raised In: Rabbit


Purity: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Storage buffer: WB: 1:500~1:1000

Form: Liquid

Reference:

Wang AH, et al. (2000) Mol Cell Biol. 20(18): 6904-6912. Grozinger CM, et al. (2000) Proc Natl Acad Sci U S A. 97(14): 7835-7840.



For research use only

Antagene, Inc. Toll Free: 1-866-964-2589, Fax: 1-888-225-1868, info@antageneinc.com,
www.antageneinc.com